

ATEC Ingenieurbüro AG
Breitenstrasse 1
6403 Küssnacht
Switzerland

Tel. +41 41 854 80 20
Fax +41 41 854 80 21
E-Mail info@swisstrac.ch
MWST-Nr. CHE-106.299.803

SWISS•TRAC®

www.swisstrac.ch

**To all
Swiss-Trac users
and transport companies**

CH-6403 Küssnacht, 18.11.2022

Official Transport Certificate and Battery Confirmation from the Manufacturer of Swiss-Trac Wheelchair Add-on Drives

Dear Swiss-Trac user,
dear transport company

The manufacturer confirms the following:

- All wheelchair tractors of the brand **Swiss-Trac** can be safely transported within any vehicle like **airplane, train, ship, public transport (e.g. bus, metro), car and others**.
The batteries must not be disassembled nor removed. No further precaution is necessary.
- The Swiss-Trac is equipped with non-spillable AGM lead batteries.
They are air- and seaworthy and comply with the IATA/ICAO special provision A67.
- The wheelchair tractor is a mobility aid that belongs to the wheelchair.
Weight mod. SWT-1/SWT-1S/SWT-1HD: approx. 70 kg / 155 lbs
Weight mod. SWT-2: approx. 55 kg / 122 lbs resp. see technical specification.

Addition for airlines:

- Airlines are obliged to carry mobility aids for people with disabilities free of charge, incl. Swiss-Trac.
- If your airline applies an additional charge, we advise you to make your reservation with another airline company instead.
- Inform your airline that the **Swiss-Trac must never be transported on the luggage belt**.
- Please use our SWISS-TRAC transport banner for air transport.

Addition for public transport companies:

- The Swiss-Trac is tip-proof.
- The wheelchair user gets more stability when the Swiss-Trac is connected to the wheelchair.

We wish you a safe and pleasant journey.

Yours sincerely,



ATEC Ing. Büro AG
Breitenstrasse 1
CH-6403 Küssnacht
Tel. +41 (0)41 854 80 20

J. Jakober, Ing. FH, Managing Director
ATEC Ingenieurbüro AG